




# Update on the ISBT TTID study on establishment of bacterial reference strains for RBC

**Dr. Marcel Prax**



**Paul-Ehrlich-Institut  
Federal Institute for Vaccines and Biomedicines  
Division Microbial Safety  
Germany**

**ISBT WP Subgroup Bacteria,  
Copenhagen, June 17<sup>th</sup>, 2017**



# Current status: overview

- ❖ **Screening for candidate strains for growth in RBC concentrates**
  - **Compilation of a panel of 7 strains for the study**
- ❖ **Production of the final bacterial solutions**
  - **Stability testing of the strains**
- ❖ **Completion of the final study protocol**
- ❖ **Coordination of participating laboratories worldwide**
- ❖ **Organization of the shipment**

# Putative bacterial candidates

## Gram positive

- *Bacillus cereus*
- *Bacillus licheniformis*
- *Listeria monocytogenes* (2x)
- *Staphylococcus epidermidis* (2x)
- *Staphylococcus aureus* (2x)
- *Streptococcus pyogenes*
- *Pseudomonas aeruginosa* (3x)
- *Pseudomonas putida*
- *Aeromonas veronii*
- *Yersinia enterocolitica* (4x)
- *Escherichia coli*
- *Morganella morganii*
- *Acinetobacter junii*
- *Klebsiella oxytoca*
- *Klebsiella pneumoniae*
- *Enterobacter cloacae*
- *Salmonella cholerasuis*

## Gram negative


- *Serratia marcescens* (3x)
- *Serratia liquefaciens* (3x)

→ 32 strains in total


# Study panel – growth behavior

strain	PEI ID	origin
<i>B. cereus</i>	PEI-B-P-57	1st WHO repository, enlargement
<i>L. monocytogenes</i>	PEI-A-199	Isolate Blood screening, England
<i>S. marcescens</i>	PEI-B-P-56	1st WHO repository, enlargement
<i>S. liquefaciens</i>	PEI-A-184	Isolate RB; C92-13-01, Roth V, et al Transfusion 2002;40(8):931-5, CDC
<i>P. fluorescens</i>	PEI-B-P-77	1st WHO repository, enlargement
<i>Y. enterocolitica</i>	PEI-A-105	Isolate RBC, Japan
<i>Y. enterocolitica</i>	PEI-A-176	Isolate RBC, CDC

*L. monocytogenes*



*S. liquefaciens*



# Stability testing

Strain	CFU/ml		
	1 week	1 month	6 month
<i>B. cereus</i>	1,23E+08	-	1,46E+08
<i>L. monocytogenes</i>	7,88E+06	8,25E+06	7,29E+06
<i>S. marcescens</i>	7,62E+06	8,11E+06	
<i>S. liquefaciens</i>	1,92E+07	2,02E+07	1,65E+07
<i>P. fluorescens</i>	1,14E+07	1,41E+07	1,10E+07
<i>Y. enterocolitica</i>	5,03E+06	6,70E+06	6,79E+06
<i>Y. enterocolitica</i>	7,33E+06	6,34E+06	6,20E+06

# Overview study design

Each of the following steps has to be performed in **triplicates**

- Shipment of **3 vials** of each strain

candidate strain



serial dilution



Inoculation  
10-25 CFU



RBC bag




Determination  
of the CFU/ml  
over time

Storage under  
routine conditions  
(2-6 °C, 42 days)


# Study design

## 1. Preparation of the final inoculum



# Study design

## 2. Verification of the final inoculum






# Study design

## 3. Spiking of red blood cell concentrate

100  $\mu$ l of last 3 dilutions in triplicates


Bacterial solution






Known concentration



10-fold serial dilution



10-25 CFU/ml




Colony counting



Storage under routine conditions  
2-6 °C, 42 days

# Study design

## 4. Sampling on days 7, 14, 21, 28, 35, 42 (?)



# Documentation

## 1. CFU counting

Bacterial strain ID	Dilution 100µl of...	Plate 1	Plate 2	Plate 3	Mean value	Bacterial count [CFU/mL]	Mean Value Vial C1 [CFU/mL]
Bag 2	Dilution 1				#DIV/0!	N/A	#DIV/0!
	Dilution 2				#DIV/0!	N/A	
	Dilution 3				#DIV/0!	N/A	
	Dilution 4				#DIV/0!	N/A	
	Dilution 5				#DIV/0!	N/A	
	Dilution 6				#DIV/0!	N/A	

## 2. Identification of bacteria

Strain:
Identification (number) of sample:
Growth after day:
Macroscopic view Colony morphology:
Microscopic view: (shape: rod, coccus)
Result of Gram-staining:
Description of identification Method (down to species level, i.e. API, PCR) (Identification panel)

## 3. Lab protocol

Test strain:			
RBC Concentrates:		Exp. Day:	
Volume			
Blood type /rh			
Control Inoculum (Dilution of stock): _____ CFU/ml (mean value)			
(CFU plate 1: _____ CFU plate 2: _____ CFU plate 3: _____)			
Result of enumeration of inoculum control: _____ CFU/ml			
Bacterial growth after storage	Sampling after	Yes (Growth)	no
	7 days		
	14 days		
	21 days		
	28 days		
	35 days		
	42 days		

## 4. Questionnaire

Please complete this questionnaire and return with the first set of completed results to allow accurate assessment.

**Study Partner: Contact**

**Name:**

Contact details:  
(Postal address, fax, phone, e-mail)

Were you a participant of the first WHO-ISBT International Validation Study on Blood Bacteria Standards? **Yes / No**

**Lab equipment used:**

Microbiological Safety Cabinet (Class II) / Laminar flow hood: **Yes / No**  
If yes, please give details: (Make, model)

# Study partners (preliminary)

Sandra Ramirez-Arcos  
Canadian Blood  
Service,  
Ottawa, Canada

Carl McDonald, Kate  
Aplin, Anjana Roy,  
NHS Blood and  
Transplant,  
London, England

Axel Seltsam, Bernd  
Lambrecht, Michael  
Schmidt, Erhard Seifried  
German Red Cross Blood  
Service NSTOB, Springe,  
Germany

Birgit Gathof,  
Melanie Störmer,  
University Hospital  
Cologne,  
Transfusion  
Medicine, Germany

Susanne Marschner,  
Shawn Keil  
Terumo BCT  
Biotechnologies,  
Lakewood, USA



Susanne Süßner,  
Austrian Red Cross,  
Blood centre Linz,  
Austria

Stephen J. Wagner,  
Holland Laboratory,  
American Red Cross,  
Blood Component  
Dep. Rockville, USA

Masahiro Satake,  
Japanese Red Cross,  
Tokyo, Japan

Julieta Rojo,  
Centro Nacional  
de la Transfusión  
Sanguínea  
Mexico

José Eduardo Levi,  
Fundação pró-  
Sangue/Hemocentro  
de São Paulo, Brasil

Ute Jentsch, Truscha Niekerk,  
Xoliswa Mpumlwana  
South African National Blood  
Service, Weltevreden Park, SA

Rita Lo, Lily Tang,  
Cheuk-kwong Lee,  
Hong Kong Red  
Cross BTS, Hong  
Kong

**Thanks to all the collaborating partners!**



# Acknowledgement



## Courtesy strains:

**Karin Janetzko**, Kirsten Fischer, DRK-BSD Mannheim

**Carl P McDonald**, NHS Blood and Transplant, London, England

**Sandra Ramirez-Arcos**, Canadian Blood Service, Ottawa

**Julieta Rojo**, Centro Nacional De La Transfusion Sanguinea, Mexico

**Masahiro Satake**, Japanese Red Cross, Tokyo

**Frank Sommer**, Universitätsklinikum Marburg

**Philipp Warnke**, Universitätsmedizin Rostock  
Universitätsklinikum, Basel

**Matthew Arduino**, Centers for Disease Control and Prevention, Atlanta

## Funding:

**German Ministry of health**

**Section microbial safety  
Paul-Ehrlich-Institut,  
Langen, Germany**

**Dr. Eva Spindler-Raffel  
Anja Schneider**

Philipp Windecker

**Julia Brachert**

Anna Marie Grothaus

PD Dr. Isabelle Bekeredjian-Ding  
(head of department)





**Thank you very much for your attention**

**Marcel Prax**

marcel.prax@pei.de